**Urea (BLT00060, BLT00061)**  
**Instrução de Uso**

REF	Código	Conteúdo
BLT00060	UREA 1000	R1: 4 x 200 mL; R2: 1 x 200 mL
BLT00061	UREA 250	R1: 4 x 50 mL; R2: 1 x 50 mL; R3: 1 x 5 mL

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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 374-7.
2. Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1838.
3. Talka, H. Schubert, G.E. Klin. Wochschr. 19; 43: 174.
4. Tiffany, T.O. Jansen, J.M. Burtis CA, Overton JB, Scott CD. Clin Chem 1972; 18: 829-40.
5. Kaplan, LA. in "Clinical Theory, Analysis and Correlation." Kaplan LA, Pesce AJ. (Ed) C V Mosby Company St Louis 1984; 1257-61.
6. Shephard, MD, Mezzachi, RD. Clin. Biochem. Revs. 1983; 4: 61-7.
7. Young, D.S. Effects of Drugs on Clinical Laboratory Tests. Third Edition. 1990; 21: 5.
8. Wachtel, M. et al, Creation and Verification of Reference Intervals. Laboratory Medicine 1995; 26: 593-7.
9. National Committee for Clinical Laboratory Standards. User evaluation of Precision Performance of Clinical Chemistry Devices. NCCLS 1984, NCCLS, Publication EP5-T.
10. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
  11. EN ISO 13485:2016 Medical devices - Quality management systems - Requirements for regulatory purposes part 2
  12. EN ISO 14971:2012 Medical devices - Application of risk management to medical devices
  13. EN ISO 15193:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for content and presentation of reference measurement procedures
  14. EN ISO 15194:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for certified reference materials and the content of supporting documentation
  15. EN ISO 15223-1:2017 Medical devices-Symbols to be used with medical devices labels, labelling and information to be supplied - Part 1: General requirements
  16. EN ISO 17511:2004 In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials

- |     |                       |                      |  |
|-----|-----------------------|----------------------|--|
| 17. | EN ISO                | 18113-1:2012         | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements       |
| 18. | EN ISO                | 18113-2:2012         | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use |
| 19. | EN ISO                | 23640:2016<br>part 2 | In vitro diagnostic medical devices - Evaluation of stability of in vitro diagnostic reagents  |
| 20. | European<br>Directive | 98/79/EC             | In Vitro Diagnostic Medical Devices Directive  |

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