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
1. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin dependent diabetes mellitus. N Engl J Med.1993; 329:977-86.
2. Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A five-year Progress Report. Clin Chem 2001; 47:1985-92.
3. Jeppson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78-89.
4. Hoelzel W, Weykamp C et al. IFCC Reference System for measurement of HbA1c in human blood and the national standardization schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:166-74.
5. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
6. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
7. EN ISO 13485:2016 part 2 Medical devices - Quality management systems - Requirements for regulatory purposes
8. EN ISO 14971:2012 Medical devices - Application of risk management to medical devices
9. EN ISO 15193:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for content and presentation of reference measurement procedures
10. EN ISO 15194:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for certified reference materials and the content of supporting documentation
11. EN ISO 15223-1:2017 Medical devices-Symbols to be used with medical devices labels, labelling and information to be supplied - Part 1: General requirements
12. EN ISO 17511:2004 In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials
13. EN ISO 18113-1:2012 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements
14. EN ISO 18113-2:2012 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use
15. EN ISO 23640:2016 part 2 In vitro diagnostic medical devices - Evaluation of stability of in vitro diagnostic reagents
16. European 2017/746 Directive In Vitro Diagnostic Medical Devices Directive

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