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REF	Código	Número de testes
REG00065	EC CARTRIDGE XS	250
REG00058	EC CARTRIDGE S	500
REG00056	EC CARTRIDGE M	1000
REG00057	EC CARTRIDGE L	3000

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1. Lewenstam, A, et al.: Application of Ion-selective Electrodes in Clinical Analysis. *Electroanalysis* 1991; 3: 727-734.
2. Lewenstam A; Routines and Challenges in Clinical application of Electrochemical Ion-Sensors. *Electroanalysis*, 2014; 26_1171-1181.
3. Cheng, S., Schneider, E. i., & Scott, M. G. (2018). Disorders of Water, Electrolytes and Acid-Base Metabolism. In N. Rifai, A. R. Horvath, & C. T. Wittwer, *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics* (pp. 1324-1347). Elsevier.
4. M. G. Scott, V.A. LeGrys, J.S. Klutts (2006). Electrolytes and blood gases. Ch. 27, in C.A. Burtis, E.R. Ashwood, & D.E. Bruns, *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 4th Ed. (pp.983-1018) Elsevier.
5. D.B. Endres, R.K. Rude, (2006) Mineral and Bone Metabolism, ch.49, in C.A. Burtis, E.R. Ashwood, & D.E. Bruns, *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 4th Ed. (pp.1891-1904) Elsevier.
6. Baird G S, Ionized Calcium. *Clin Chim Acta* 2011; 412: 696-701.
7. Burnett RW, Christiansen TF, Covington AK, et al. IFCC recommended reference method for the determination of the substance concentration of ionized calcium in undiluted serum, plasma, or whole blood. *Clin Chem Lab Med*. 2000; 38 (12): 1301-1314. doi: 10.1515/CCLM.2001.206
8. IFCC 2000/1; Burnett RW, Covington AK, Fogh-Andersen N, et al., Use of Ion-Selective Electrodes for Blood-Electrolyte Analysis. Recommendations for Nomenclature, Definitions and Conventions, *Clin Chem Lab Med* 2000; 38 (4): 363-370.
9. IFCC 2000/3 Burnett RW, Covington AK, Fogh-Andersen N, et al. Recommendations for measurement of and conventions for reporting sodium and potassium by Ion-selective electrodes in undiluted serum, plasma, or whole blood. *Clin Chem Lab Med*. 2000; 38 (10): 1065-1071. doi: 10.1515/CCLM.2000.159
10. Ben Rayana MC, Burnett RW, Covington AK, et al. IFCC Recommendations for measuring and reporting chloride by ISEs in undiluted serum, plasma, or blood. *Clin Chem Lab Med*. 2006; 44 (3): 346-352. doi: 10.1515/CCLM.2006.060
11. CLSI. *Ionized Calcium Determinations: Precollection Variables, Specimen Choice, Collection, and Handling; Approved Guideline - Second Edition*. CLSI document C31-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2001.
12. Boink, A B T J, et al, IFCC recommendations on sampling, transport, and storage for the determination of the concentration of ionized calcium in whole blood, plasma, and serum. *J Autom. Chem*, 1991; 13 (5): 235-239.
13. P. Buhlmann and L. D. Chen, "Ion-Selective Electrodes with Ionophore-Doped Sensing Membranes", in *Supramolecular Chemistry*, P.A. Gale and J. W. Steed, Eds. 2012, pp. 2539-2579.
14. E. Lindner and B. D. Pendley, A tutorial on the application of ion-selective electrode potentiometry: An analytical method with unique qualities, unexplored opportunities, and potential pitfalls. *Anal. Chim. Acta*, vol. 762, pp. 1-13, 2013.
- 15 EN ISO 13485:2016 part 2 Medical devices - Quality management systems - Requirements for regulatory purposes
- 16 EN ISO 14971:2012 Medical devices - Application of risk management to medical devices
- 17 EN ISO 15193:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for content and presentation of reference measurement procedures
- 18 EN ISO 15194:2009 In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for certified reference materials and the content of supporting documentation
- 19 EN ISO 15223-1:2017 Medical devices-Symbols to be used with medical devices labels, labelling and information to be supplied - Part 1: General requirements
- 20 EN ISO 17511:2004 In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials
- 21 EN ISO 18113-1:2012 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements
- 22 EN ISO 18113-2:2012 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use
- 23 EN ISO 23640:2016 part 2 In vitro diagnostic medical devices - Evaluation of stability of in vitro diagnostic reagents
- 24 European 98/79/EC In Vitro Diagnostic Medical Devices Directive

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