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| REF | Código | Número de testes |
|----------|--------------------------|------------------|
| REG00066 | EC Cartridge XS Plus iCa | 250 |
| REG00061 | EC Cartridge S Plus iCa | 500 |
| REG00062 | EC Cartridge M Plus iCa | 1000 |
| REG00063 | EC Cartridge L Plus iCa | 3000 |

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- Lewenstam, A, et al.: Application of Ion-selective Electrodes in Clinical Analysis. *Electroanalysis* 1991; 3: 727-734.
- Lewenstam A; Routines and Challenges in Clinical application of Electrochemical Ion-Sensors. *Electroanalysis*, 2014; 26_1171-1181.
- Cheng, S., Schneider, E. i., & Scott, M. G. (2018). Disorders of Water, Electrolytes and Acid-Base Metabolism. In N. Rifai, A. R. Horvath, & C. T. Wittwer, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics (pp. 1324-1347). Elsevier.
- M. G. Scott, V.A. LeGrys, J.S. Klutts (2006). Electrolytes and blood gases. Ch. 27, in C.A. Burtis, E.R. Ashwood, & D.E. Bruns, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th Ed. (pp.983-1018) Elsevier.
- D.B. Endres, R.K. Rude, (2006) Mineral and Bone Metabolism, ch.49, in C.A. Burtis, E.R. Ashwood, & D.E. Bruns, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th Ed. (pp.1891-1904) Elsevier.
- Baird G S, Ionized Calcium. *Clin Chim Acta* 2011; 412: 696-701.
- Burnett RW, Christiansen TF, Covington AK, et al. IFCC recommended reference method for the determination of the substance concentration of ionized calcium in undiluted serum, plasma, or whole blood. *Clin Chem Lab Med.* 2000; 38 (12): 1301-1314. doi: 10.1515/CCLM.2001.206
- IFCC 2000/1; Burnett RW, Covington AK, Fogh-Andersen N, et al., Use of Ion-Selective Electrodes for Blood-Electrolyte Analysis. Recommendations for Nomenclature, Definitions and Conventions, *Clin Chem Lab Med* 2000; 38 (4): 363-370.
- IFCC 2000/3 Burnett RW, Covington AK, Fogh-Andersen N, et al. Recommendations for measurement of and conventions for reporting sodium and potassium by Ion-selective electrodes in undiluted serum, plasma, or whole blood. *Clin Chem Lab Med.* 2000; 38 (10): 1065-1071. doi: 10.1515/CCLM.2000.159
- Ben Rayana MC, Burnett RW, Covington AK, et al. IFCC Recommendations for measuring and reporting chloride by ISEs in undiluted serum, plasma, or blood. *Clin Chem Lab Med.* 2006; 44 (3): 346-352. doi: 10.1515/CCLM.2006.060
- CLSI. *Ionized Calcium Determinations: Precollection Variables, Specimen Choice, Collection, and Handling; Approved Guideline - Second Edition.* CLSI document C31-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2001.
- Boink, A B T J, et al, IFCC recommendations on sampling, transport, and storage for the determination of the concentration of ionized calcium in whole blood, plasma, and serum. *J Autom. Chem,* 1991; 13 (5): 235-239.
- P. Buhlmann and L. D. Chen, "Ion-Selective Electrodes with Ionophore-Doped Sensing Membranes", in *Supramolecular Chemistry*, P.A. Gale and J. W. Steed, Eds. 2012, pp. 2539-2579.
- E. Lindner and B. D. Pendley, A tutorial on the application of ion-selective electrode potentiometry: An analytical method with unique qualities, unexplored opportunities, and potential pitfalls. *Anal. Chim. Acta*, vol. 762, pp. 1-13, 2013.
- EN ISO 13485:2016 part 2 Medical devices - Quality management systems - Requirements for regulatory purposes

| | | | |
|----|--------------------|-------------------|--|
| 16 | EN ISO | 14971:2012 | Medical devices - Application of risk management to medical devices |
| 17 | EN ISO | 15193:2009 | In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for content and presentation of reference measurement procedures |
| 18 | EN ISO | 15194:2009 | In vitro diagnostic medical devices - Measurement of quantities in samples of biological origin - Requirements for certified reference materials and the content of supporting documentation |
| 19 | EN ISO | 15223-1:2017 | Medical devices-Symbols to be used with medical devices labels, labelling and information to be supplied - Part 1: General requirements |
| 20 | EN ISO | 17511:2004 | In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials |
| 21 | EN ISO | 18113-1:2012 | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements |
| 22 | EN ISO | 18113-2:2012 | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use |
| 23 | EN ISO | 23640:2016 part 2 | In vitro diagnostic medical devices - Evaluation of stability of in vitro diagnostic reagents |
| 24 | European Directive | 98/79/EC | In Vitro Diagnostic Medical Devices Directive |

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